

7 LED COLORS FACE THERAPY WAND USER MANUAL



zeencare.com support@zeencare.com

INCLUSIONS

- Face Wand
- 2 Charging Cable
- 3 User Manual

PRODUCT SPECS

- Rotating Head
- Power Button/ Color Button
- EMS Intensity Level Button
- EMS LED Indicator



HOW TO USE

- Power On/Off: Press and hold the upper button to turn on/off the Face Wand.
- Color Selection: Quickly press the upper button to change colors.
- EMS Intensity: Press and hold the lower button to turn on/off EMS. Press again to adjust intensity from low to high.
- Use on Clean Skin: Ensure your skin is clean and moisturized before using.
- Care and Storage: Wipe the wand after use and store it in a dry place away from sunlight.

7 LED COLORS







Blue

Cyan

White

Green

Yellow

Purple

PAGE 4

PRODUCT SPECS

Face Wand Size	15L x 2.8W x 2.2H 38cmL x 7cmW x 5.5 cmH
Product Weight	43g/0.095 lbs
Material	Aluminum Alloy
Charging Type	Туре С
Charging Time	3hrs

- Warm Therapy
- Mild Microcurrent
- Ergonomic Shape Head Massage
- 3 EMS Instensity Levels
- LED Light Indicator
- Auto Shutdown after 10 mins

Tips

- Charge the battery fully for the first use.
- Clean equipment with a dry cloth to remove residue and keep it dry.
- Use moisturizer or any of your preferred face cream when using the Face Wand.
- Don't press too hard on the skin.
- Adjust intensity gradually to your comfort level.
- Use the wand regularly for best results.
- Clean the wand after each use.
- Store it in a dry place away from direct sunlight.

- Don't use on damaged or irritated skin.
- Prevent direct contact with eyes.
- Avoid water exposure to prevent damage.
- Don't use product while charging.
- Don't disassemble the Face Wand

Enhance Your Beauty Routine

Welcome to the world of face wands – your gateway to radiant, youthful skin!

From reducing puffiness and improving circulation to boosting collagen production and minimizing fine lines, the Zeencare face wand is here for every skin concern and goal.

Revolutionize your skincare routine, leaving you with glowing, rejuvenated skin you'll love to flaunt.