

SPECIFICATION

Product Name	Under Eye Red Light Therapy
Product Model	ZC-LIT-110670
Features	Microcurrent, Red Light Therapy, 4 Therapy Modes, 10-min Timer, Rechargeable Battery
Battery Capacity	90 mAh
Rated Voltage	3.7 V
Included Components	Under Eye Therapy Device, USB Type-C Cord, Manual
Product Dimension	5.5 x 5.9 x 1.7 in 14 x 15 x 4.5 cm
Product Weight	1.34 oz 38 g

DESCRIPTION

Combines gentle Microcurrent and Red Light Therapy to help reduce under-eye puffiness, dark circles, and fine lines. Designed for daily use, it revitalizes the under-eye area for a smoother, brighter, and refreshed appearance.

Zeencare UNDER EYE RED LIGHT THERAPY



Zeencare

Promoting Health, Inspiring Wellness

✉ support@zeencare.com

🌐 www.zeencare.com

USAGE GUIDE



- Long Press: Power On/Off
Short Press: Mode Switch
- Microcurrent Intensity
Light Indicator
- Adjust Microcurrent Intensity
Short Press: Low, Medium, High

MODE SWITCH INSTRUCTIONS

1st Short Press – Mode 1: Revitalize & Firm

Microcurrent fixed frequency + Red Light Therapy

2nd Short Press – Mode 2: Lift & Illuminate

Microcurrent variable frequency + Red Light Therapy

3rd Short Press – Mode 3: Tone & Relax

Microcurrent variable frequency

4th Short Press – Mode 4: Refresh & Renew

Red Light Therapy only

SAFETY NOTES & USAGE TIPS

1. Ensure your face is clean and free of makeup before using.
2. Apply a skin gel or moisturizer for better results.
3. A slight tingling sensation is normal with Microcurrent devices. Sensation varies with skin moisture – more moisture may increase intensity.
4. This device is not waterproof. Wipe with a damp cloth; do not rinse with water.

CHARGING GUIDE

- Charging: Red light stays on while charging. Green light indicates fully charged.
- Low Battery Reminder: Red light flashes 3 times when the battery is low.

